

GSPM Table™

www.medicalpositioning.com
800-593-3246



 **MPI**
Positioned for results

STANDARD FEATURES

| | | | | |
|------------------------------|-----------------------|-----------------------|----------------------|------------------------|
| Height Adjustability: | Load Capacity: | Lift Capacity: | Imaging Area: | Attenuation: |
| 29" to 39" | 550 lbs | 550 lbs | 53" | 0.57 mm Al. equivalent |

| | |
|------------------------|--|
| Control System: | Additional Features: |
| Hand Controller | Cervical Access Shaped Top Prone Procedure Facial Opening Patient Restraints C-Arm Accessible Base Carbon Fiber Top and 2" Pad |

OPTIONS

| |
|--|
| Trendelenburg: +20° to -20° |
| Lateral Roll: -20° to +20° |
| Hand Controller with Memory Positioning & Battery Backup |
| Footswitch |
| Extended Warranty |

PRODUCT SPECIFICATIONS

| | | |
|----------------|---------------|----------------|
| Length: | Width: | Weight: |
| 84" | 22" | 290 lbs |

| | | |
|-------------------|---------------------------------|---|
| Foam: | Electrical: | Vinyl: |
| Cal. B.F.T.B #117 | 120 VAC, 7.0 amps max, 50/60 Hz | Fed. Spec. Cec-A-680A D.O.T. FAR 25.8536, M.V.S. 302 Port of NY/Boston F.D. Code |

UL for Hospital Use:

UL/IEC 60601-1, CAN/CSA C22.2 No. 601.1

Warranty:

5 Years - All electrical, mechanical, and structural; 1 Year parts and labor (see Warranty for complete details)

ACCESSORIES

| | |
|---|---|
| Product Specific: | Universal: |
| Dual Articulating Arm Boards Accessory Rail Extender | 816 or 624 Head Elevation Safe-T-Wedge 820 or 626 Body Positioning Safe-T-Wedge Set of 2 Safety Rails Paper Roll Holder & Cutter |

COLOR OPTIONS

| | |
|------------------|---|
| Standard: | Available with additional lead time: |
| ● Space Blue | ● Spanish Red ● Desert Sage ● Taupe ● Oatmeal |