


HUT Table™

www.medicalpositioning.com
800-593-3246



Positioned for results

STANDARD FEATURES

Height Adjustability:	Load Capacity:	Tilt Rate:	Trendelenburg:
30" to 38"	500 lbs	6.5°/second	+85° to -15°
Control System:	Additional Features:		
Hand Controller with Battery BackUp	Dual Fail Safe Back Up System Built-in Folding Foot Platform Patient Restraints		

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
73"	30"	380 lbs
Foam:	Electrical:	Vinyl:
Cal. B.F.T.B #117	120 VAC, 1.6 amps max, 50/60 Hz	D.O.T. FAR 25.8536, M.V.S. 302 Port of NY/Boston F.D. Code

Industry Standards:

UL/IEC 60601-1, CAN/CSA C22.2 No. 601.1

Warranty:

1 Year - All electrical, mechanical, and structural; 1 Year parts and labor (see Warranty for complete details)






OPTIONS

Drop Sections - Imaging Drop Section - Sonographer Drop Section	Footswitch
---	------------

ACCESSORIES

Product Specific:	Universal:
Pediatric/Geriatric Adapter HUT Head Extension Fixed Arm Boards L Shaped Arm Boards Carotid Head Support	816 or 624 Head Elevation Safe-T-Wedge 820 or 626 Body Positioning Safe-T-Wedge Set of 2 Safety Rails

COLOR OPTIONS

Standard:	Available with additional lead time:
 Bay Leaf	 Spanish Red  Desert Sage  Taupe  Oatmeal