



Positioned for results

STANDARD FEATURES

Height Adjustability:	Load Capacity:	Drop Sections:	Trendelenburg:	Braking:
24" to 34"	1,000 lbs	Imaging Drop Section Sonographer Drop Section	+15° to -15°	Single Pedal Braking

Control System:	Ergonomical Features:	Additional Features:	OPTIONS
Hand Controller with Battery BackUp	Sonographer Extension 2-Way Sonographer Drop Section Patient Back Support	Safety Rails Paper Roll Holder & Cutter	65° Fowler (Head Up) Positioning
			Always EchoBed® (No Fowler)
			Footswitch
			Upgrades to Stress EchoBed® purchase

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
74"	30"	340 lbs
Foam:	Electrical:	Vinyl:
Cal. B.F.T.B #117	120 VAC, 1.6 amps max, 50/60 Hz	D.O.T. FAR 25.8536, M.V.S. 302 Port of NY/Boston F.D. Code

Industry Standards:

UL/IEC 60601-1, CAN/CSA C22.2 No. 601.1

Warranty:

5 Years - All electrical, mechanical, and structural; 1 Year parts and labor (see Warranty for complete details)

ACCESSORIES

Product Specific:	Universal:
Pediatric/Geriatric Adapter Padded Arm Boards Leg Supports IV Pole + Holder Carotid Head Support	816 or 624 Head Elevation Safe-T-Wedge 820 or 626 Body Positioning Safe-T-Wedge

COLOR OPTIONS

Standard:	Available with additional lead time:
● Bay Leaf	● Spanish Red ● Desert Sage ● Taupe ● Oatmeal