


Stress EchoBed®

www.medicalpositioning.com
800-593-3246



Positioned for results

STANDARD FEATURES

Height Adjustability:	Load Capacity:	Drop Sections:	Trendelenburg:	Braking:
24" to 34"	1,000 lbs	Imaging Drop Section Sonographer Drop Section	+15° to -15°	Single Pedal Braking
Control System:	Ergonomical Features:		Additional Features:	
Hand Controller with Memory Positioning & Battery Backup	Sonographer Extension 2-Way Sonographer Drop Section Patient Back Support		Safety Rails Lateral Tilt: 0° to 40°	Ergometer with display designed for stress echo including 10 preset exercise echo protocols, 1 manual protocol and 4 programmable protocols

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
78"	30"	465 lbs
Foam:	Electrical:	Vinyl:
Cal. B.F.T.B #117	120 VAC, 3.2 amps max, 50/60 Hz	D.O.T. FAR 25.8536, M.V.S. 302 Port of NY/Boston F.D. Code

Industry Standards:

UL/IEC 60601-1, CAN/CSA C22.2 No. 601.1

Warranty:

5 Years - All electrical, mechanical, and structural; 1 Year parts and labor
(see Warranty for complete details)

ACCESSORIES

Product Specific:

Pediatric/Geriatric Adapter
Padded Arm Boards
IV Pole + Holder
Carotid Head Support
Storage Cart & Track Pad
ECG Workstation Interface
Ergometer Controller
Ergometer Insert

Universal:

816 or 624 Head Elevation
Safe-T-Wedge
820 or 626 Body Positioning
Safe-T-Wedge

COLOR OPTIONS

Standard:

● Bay Leaf

Available with additional lead time:

● Spanish Red

● Desert Sage

● Taupe

● Oatmeal