

UltraScan Table™

www.medicalpositioning.com
800-593-3246



 **MPI**
Positioned for results

STANDARD FEATURES

Height Adjustability:	Load Capacity:	Fowler (Head Up):	Foot Down Section:	Braking:
24.5" to 34.5"	500 lbs	0° to 70°	Manual Adjustment	Individual Locking Casters
Control System:	Additional Features:			
Hand Controller	Removable Foot Drop Section Adjustable Stirrups			

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
72"	28"	290 lbs
Foam:	Electrical:	Vinyl:
Cal. B.F.T.B #117	120 VAC, 1.6 amps max, 50/60 Hz	D.O.T. FAR 25.8536, M.V.S. 302 Port of NY/Boston F.D. Code

Industry Standards:

UL/IEC 60601-1, CAN/CSA C22.2 No. 601.1

Warranty:

1 Year - All electrical, mechanical, and structural; 1 Year parts and labor (see Warranty for complete details)

OPTIONS

Dual Drop Sections - Imaging Drop Section - Sonographer Drop Section	Hand Controller with Memory Positioning & Battery Backup
Trendelenburg +15° to - 15°	Footswitch
Trendelenburg +25° to - 15°	Fluid Containment
Single Pedal Braking	

ACCESSORIES

Product Specific:	Universal:
Pediatric/Geriatric Adapter Padded Arm Boards Leg Supports IV Pole + Holder Carotid Head Support	816 or 624 Head Elevation Safe-T-Wedge 820 or 626 Body Positioning Safe-T-Wedge Set of 2 Safety Rails Paper Roll Holder & Cutter

COLOR OPTIONS

Standard:	Available with additional lead time:
● Bay Leaf	● Spanish Red ● Desert Sage ● Taupe ● Oatmeal