

# VasScan Table™

[www.medicalpositioning.com](http://www.medicalpositioning.com)  
800-593-3246



 **MPI**  
Positioned for results

## STANDARD FEATURES

<b>Height Adjustability:</b>	<b>Load Capacity:</b>	<b>Fowler (Head Up):</b>	<b>Trendelenburg:</b>	<b>Braking:</b>
26.5" to 34.5"	500 lbs	0° to 60°	+30° to -15°	Individual Locking Casters
<b>Control System:</b>	<b>Additional Features:</b>			
Hand Controller	Adjustable and Removable Foot Platform			

## PRODUCT SPECIFICATIONS

<b>Length:</b>	<b>Width:</b>	<b>Weight:</b>
80"	28"	265 lbs
<b>Foam:</b>	<b>Electrical:</b>	<b>Vinyl:</b>
Cal. B.F.T.B #117	120 VAC, 1.6 amps max, 50/60 Hz	D.O.T. FAR 25.8536, M.V.S. 302 Port of NY/Boston F.D. Code

### Industry Standards:

UL/IEC 60601-1, CAN/CSA C22.2 No. 601.1

### Warranty:

1 Year - All electrical, mechanical, and structural; 1 Year parts and labor (see Warranty for complete details)

## OPTIONS

Dual Drop Sections - Imaging Drop Section - Sonographer Drop Section	Footswitch
Height Adjustability: 28.5" - 36.5"	Hand Controller with Memory Positioning & Battery Backup
Height Adjustability: 30.5" - 38.5"	
Single Pedal Braking	

## ACCESSORIES

<b>Product Specific:</b>	<b>Universal:</b>
Pediatric /Geriatric Adapter	816 or 624 Head Elevation Safe-T-Wedge
Padded Arm Boards	820 or 626 Body Positioning Safe-T-Wedge
Leg Supports	Set of 2 Safety Rails
IV Pole & Holder	Paper Roll Holder & Cutter
Carotid Head Support	

## COLOR OPTIONS

<b>Standard:</b>	<b>Available with additional lead time:</b>
<span style="color: #8ebf42;">●</span> Bay Leaf	<span style="color: #c00000;">●</span> Spanish Red <span style="color: #556b2f;">●</span> Desert Sage <span style="color: #808080;">●</span> Taupe <span style="color: #f5deb3;">●</span> Oatmeal