


— DBI Table™

The ultimate platform for Decubitus Stereotactic Breast Biopsy

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246



ergonomic * innovative * proven

Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Drop Sections:	Control System:	Braking:
26" to 38"	500 lbs	Four imaging drop sections; two on each side to allow space for the image receptor and physician	Hand Controller	Individual Locking Casters

Additional Features:

Back Support
Interchangeable Head and Foot Sections

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
75"	27"	315 lbs
Foam:	Electrical:	Vinyl:
California Technical Bulletin 117 - Section E	120 VAC, 50/60 Hz, 2.0A Max 230 VAC, 50/60 Hz, 1.0A Max	FMVSS 302 Boston Fire Code - BFD IX-I NY/NJ Port Authority

CERTIFICATIONS

Industry Standards:	Warranty:
UL 60601-1, CAN/CSA C22.2 No. 601.1	1 Year

PRODUCT BENEFITS

360° access to breast
Increased access to lesions in upper/outer quadrant
Compatible with any stereotactic add-on system
Drop sections provide closer access for image receptor
Decubitus positioning provides greater patient comfort
Enhances ergonomic posture as technicians work upright in drop sections
Back support allows the patient to take deep breaths while keeping the chest wall in contact with the image receptor
Minimizes risks associated with fainting during an upright procedure
Bariatric capable

OPTIONAL FEATURES

Foot Controller