


# GSPM Table™

The ultimate platform for General Surgery and Pain Management

MEDICAL POSITIONING, INC.  
www.medicalpositioning.com  
800-593-3246



ergonomic \* innovative \* proven

Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

**STANDARD FEATURES**

<b>Height Adjustability:</b>	<b>Maximum Patient Capacity:</b>	<b>Attenuation:</b>	<b>Trendelenburg/Reverse Trendelenburg:</b>	<b>Lateral Tilt:</b>
29" to 39"	550 lbs	0.57 mm Aluminum Equivalent	20°/20°	20° Left and Right
<b>Control System:</b>	<b>Braking:</b>	<b>Imaging Area:</b>	<b>Additional Features:</b>	
Hand Controller	Individual Locking Casters	53"	Cervical Access Shaped Top Prone Procedure Facial Opening C-Arm Accessible Base	
			Carbon Fiber Top and 2" Pad Patient Safety Straps	

**PRODUCT SPECIFICATIONS**

<b>Length:</b>	<b>Width:</b>	<b>Weight:</b>
84"	22"	290 lbs
<b>Foam:</b>	<b>Electrical:</b>	<b>Vinyl:</b>
California Technical Bulletin 117 - Section E	120 VAC, 60 Hz, 1.6A Max	FMVSS 302 Boston Fire Code - BFD IX-I NY/NJ Port Authority

**CERTIFICATIONS**

<b>Industry Standards:</b>	<b>Warranty:</b>
UL 60601-1, CAN/CSA C22.2 No. 601.1	5 Years

**PRODUCT BENEFITS**

Low attenuation reduces scatter for better imaging  
 Cantilever design for better C-Arm accessibility during pain management procedures  
 Bariatric capable

**OPTIONAL FEATURES**

Foot Controller  
 Carbon Fiber Arm Boards