



Convertible platform for 360° access to breast

## Chair to table positioning for multi use

The UltraDBI maximizes department space by incorporating upright and decubitus stereotactic breast biopsy, as well as mammography into one bariatric capable platform, increasing throughput and increasing continuity of care.

- Easy transition from chair to table provides an all-encompassing platform for varying patient population
- Accommodates both upright and decubitus positions for stereotactic breast biopsy, as well as mammography and breast ultrasound
- Bariatric capable patient weight capacity of 500 lbs establishes continuity of care
- Reduce turnover time and increase throughput, even in small exam rooms, with interchangeable head and foot rests eliminating the need to reposition the table to accommodate each patient
- Accessories Include: removable arm rests, back support, interchangeable, adjustable and removable head supports

**STANDARD FEATURES**

<b>Height Adjustability:</b>	<b>Maximum Patient Capacity:</b>	<b>Trendelenburg:</b>	<b>Fowler:</b>	<b>Drop Sections:</b>	<b>Supported Patient Positions:</b>
25" to 34.5"	500 lbs	15° to -15°	90°	1 on each side: one for physician and one for receptor	Upright in a chair and Decubitus on a table
<b>Foot Platform:</b>	<b>Control System:</b>	<b>Braking:</b>	<b>Additional Features:</b>		
Adjustable foot platform in the chair position	Hand Controller	Single pedal braking	Adjustable Arm Rest	Detachable Head Support	Adjustable Back Support Adjustable Head Support

**PRODUCT SPECIFICATIONS**

<b>Length:</b>	<b>Surface Width:</b>	<b>Length of Wheel Base:</b>	<b>Foam:</b>	<b>Electrical:</b>	<b>Vinyl:</b>
76"	28"	30"	California Technical Bulletin 117 - Section E	100-230 VAC, 50/60 Hz, Max 5 A	FMVSS 302 Boston Fire Code - BFD IX-I NY/NJ Port Authority

<b>Certifications:</b>	<b>Warranty:</b>	<b>OPTIONAL FEATURES</b>			
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	1 Year	IV Pole and Holder Battery	Foot Control	Back Support Wedge Cushion	Underbed Lights
		Instrument Tray			